
L&L SKIN

MIKO
User Manual

Discover radiant with the MIKO, which features EP, heating, and red,blue, purple, and yellow light therapy, all complemented by gentle vibrations. This innovative device provides a luxurious and heavenly experience all in one. The advanced technology will enhance your complexion's radiance while boosting the absorption rate of serums and moisturizers.Experience the transformative effects of MIKO and elevate your self-care routine to new heights



PRODUCT MODE:

Glow EP+Red LED light

Enhance skin radiance and collagen production with EP currents and red light for a youthful glow.

Sculpt EP+Blue LED light

Define facial contours and rejuvenate skin by tightening and reducing dark spots.

Restore heating+Purple LED light

Boost blood circulation, combat inflammation, and improve skin elasticity for a flawless complexion.

Massage Vibration+Yellow LED light

Relax and rejuvenate skin, reducing fatigue and oxidation to slow aging

How to use



Watering
Place the head against the skin and roll it in a circular motion, providing a boosting effect.



Boosting
Press and slide the head against the skin as if gliding, providing a boosting effect.



Focus
Place the head against the skin and press it on one spot for 1-2 seconds before moving to a next spot

Cleaning suggestions:
*Clean immediately after each use to avoid dirt and bacteria accumulation,
*Avoid soaking. Do not immerse the beauty instrument in water or place it in a humid environment.

Japan L&L Skin Tech Co., Limited
3F Kabutocho Daiichi Heiwa Building5-1
Nihonbashi Kabuto-cho, Chuo-ku, Tokyo, Japan Post code:103-026

www.llskin.jp
www.llskinshop.com

EAC CE RoHS ISO